

Efficient work processes and easy handling

ZEISS SEM Dual Joystick and Control Panel



Sample picture of a Control Panel.

Upgrade your ZEISS microscope

Upgrade your ZEISS microscope with the SEM Dual Joystick and Control Panel and benefit from easy handling, more comfort and efficient work processes.

As technical requirements may apply on some systems, please contact us to learn more about the SEM Dual Joystick and Control Panel and how your processes will benefit from an upgrade:
microscopy@zeiss.com

There are several ways you can operate Scanning Electron Microscopes (SEM). With the Dual Joystick Controller and the Control Panel installed, operating becomes much more comfortable. The Dual Joystick Controller can be used for stage control and specimen navigation, and the Control Panel gives you easy access to the most frequently used functions of the SEM.

Highlights

Various configuration options

Use the high-resolution encoders on the full-sized keyboard of the Control Panel to control functions such as focus, magnification, contrast and brightness

Many compensation possibilities

The Dual Joystick allows sensitive and manual positioning.
All axes are deflection- and magnification-compensated

Convenient handling

Save time thanks to free and precise adjustment of travel speed and easy cleaning



Sample picture of a Dual Joystick.



Seeing beyond

ZEISS SEM Dual Joystick and Control Panel

Easy handling



Availability

- Crossbeam series
- GeminiSEM series
- MERLIN series
- Sigma series
- EVO series
- AURIGA series
- ULTRA series
- SUPRA series
- NVision series
- NEON series
- 15xx XB series
- 15xx series



microscopy@zeiss.com
www.zeiss.com/microscopy