

## Optical Biometry Educational Event

<b>09.30 – 10.00</b>	<b>Coffee and Registration</b>
<b>10.00 – 11.00</b>	<b>Introduction</b> <ul style="list-style-type: none"><li>• Recap Optical Biometry</li><li>• Refraction refresher</li><li>• Tips &amp; Tricks</li><li>• RCO Guidelines</li></ul>
<b>11.00 – 11.15</b>	<b>Coffee</b>
<b>11.15 -12.45</b>	<b>Workshop Rotation</b> <ul style="list-style-type: none"><li>• IOL Master workflow</li><li>• Tips &amp; Tricks for challenging situations</li><li>• IOL Master 700</li></ul>
<b>12.45- 13.30</b>	<b>Lunch</b>
<b>13.30 - 14.00</b>	Formulae and Lens Constant /Optimisation
<b>14.00 – 15.30</b>	<b>Workshops</b> <ul style="list-style-type: none"><li>• Advanced Tips and Tricks</li><li>• Case Studies</li><li>• Loading Surgeons, Lens Constants/Optimisation screen</li></ul>
<b>15.30 – 15.50</b>	<b>Presentations</b> <ul style="list-style-type: none"><li>• Complete Cataract Solution</li><li>• IOL Master 700</li><li>• Feedback &amp; Validation</li></ul>
<b>16.00</b>	<b>Close</b>