

Sl.No	Segment	Product
1	Surgical Micorscopes	OPMI LUMERA 300
		OPMI LUMERA i
		OPMI LUMERA 700
		ARTEVO 800
		Callisto eye
2	Biometry	IOLMaster 500
		IOLMaster 700
3	Phacoemulsification	VISALIS 100
		VISALIS 500
4	IOLs	IOL portfolio
5	miLOOP	miLOOP
6	FCI	FCI
7	Slit Lamps	SL 115 / 120 / 130 / 220
		SL 800
		SL Imaging
8	Hand held Fundus camera	VISUSCOUT 100
9	Autorefractor	VISUREF 150
10	Lensmeter	VISULENS 500/550
11	OCT	CIRRUS 500 / 5000 Angioplex
		CIRRUS 6000
12	Perimeter	HFA 3
13	Fundus imaging	CLARUS 500 / 700
		VISUCAM 524
14	Ophthalmic laser	VISULAS green
		VISULAS YAG III
15	Connectivity	FORUM
16	Refractive	MEL 90
		VISUMAX
		SMILE
		PRESBYOND