ZEISS Medical Technology Open House

Quality Solutions for Medical Standards

Wednesday, May 25th ZEISS Quality Excellence Center Wixom, MI



Manage your medical regulatory requirements

Looking for ways to stay ahead or need a competitive edge in the ever-changing medical industry?

Expand your quality solutions knowledge with ZEISS experts in person! You'll see demonstrations of the latest ZEISS technologies and hear from those who work with your measurement applications every day.

Learn how to:

- Solve new inspection challenges brought on by medical requirements and new manufacturing methods
- Apply the proper inspection technologies to optimize your measurement processes
- Leverage quality data to gain real-time insights and drive process improvements

Register Today

Agenda:

8:00 am - 9:00 am Welcome, Pastries, Coffee 9:00 am - 3:00 pm Visit technology stations and demos:

- CT X-ray Technologies for Medical Devices
- Tactile, Optical, and White Light Applications
- 3D Scanning and Inspection with GOM Technology
- Compliant Microscopy Processes
- Software Validation and FDA Compliance
- Technical Cleanliness
- Additive Manufacturing in Medical Applications

12:00 pm - 1:00 pm Lunch provided by ZEISS

Join us at any point in the day that works for you! Please register to save your spot by May 18th.

Event Address:

ZEISS Quality Excellence Center 29295 Lyon Oaks Dr. Wixom, MI 48393



Attendees will be entered for a chance to win a Bose SoundLink portable speaker!



Seeing beyond