

Segment	Shop Online	Apply for EMI offer
Surgical	ZEISS OPMI LUMERA 300 Basic System	-
	ZEISS IOLMaster 700	-
	CT & TK Upgrade for ZEISS IOLMaster 700	Apply here
	ZEISS miLOOP	-
Diagnostics & Therapeutics	HFA3 840	-
	CIRRUS HD-OCT 500	-
	CIRRUS HD OCT 500 (Demo Unit)	-
Data Management	ZEISS FORUM Office	-
	ZEISS Glaucoma Workplace	
	ZEISS EQ Workplace	
	ZEISS Retina Workplace	
Essential line	VISULENS 550 - Digital Lensmeter	
	VISUREF 150 - Autorefractor/Keratometer	
	VISUPLAN 500 - Non-contact Tonometer	
	VISUSCOUT 100 - Handheld Fundus Camera	Apply here
	ZEISS SL 115 Slit Lamp	Apply here
	Applanation Tonometer AT 020 For SL 115 Classic	
	ZEISS SL 220 / 3-Step Magnification	
	Applanation Tonometer AT 030 - Swiveling for SL 220	
FCI Vitreoretinal products	Surgical Retina Medical Retina	
Eye Centration Device	i.Terminal go	
Laboratory microscope	ZEISS Primostar 1	
	ZEISS Primostar 3	
	ZEISS Stemi 508	